



# PRODUCT LIST 2021

Immunoassays ELISA



## REPRODUCTIVE

Kit	Cat #	Wells	Enzyme	Wave-length
AMH Gen II	A79765	12x8	HRP	450
AMH Gen II Cal+CS	A79766	N.A.	N.A.	N.A.
AMH Gen II RUO	A73818	12x8	HRP	450
AMH Gen II Cal+CS RUO	A73819	N.A.	N.A.	N.A.
Inhibin A	DSL-10-28100t-1	12x8	HRP	450
Inhibin A	DSL-10-28100t-4	48x8	HRP	450
Inhibin B Gen II	A81303	12x8	HRP	450
Inhibin B Gen II Cal+CS	A81304	N.A.	N.A.	N.A.
Inhibin B Gen II RUO	A81301	12x8	HRP	450
Inhibin B Gen II Cal+CS RUO	A81302	N.A.	N.A.	N.A.
Free Testosterone	<b>NEW</b> C52296	12x8	HRP	450

Kit	Sensitivity (LoD)	Incubation time	Calibrator range	Sample	
				Vol. [µL]	Form
AMH Gen II	0.05 ng/mL	2x1 hr/RT; SH+30+10 min/RT; SH	0; 0.16-22.5 ng/mL	30	s; p
AMH Gen II Cal+CS	N.A.	N.A.	N.A.	N.A.	N.A.
AMH Gen II RUO	0.05 ng/mL	2x1 hr/RT; SH+30+10 min/RT; SH	0; 0.16-22.5 ng/mL	30	s; p
AMH Gen II Cal+CS RUO	N.A.	N.A.	N.A.	N.A.	N.A.
Inhibin A	2.61 pg/mL	3 hrs+1 hr+15 min/RT; SH	0; 10-1000 pg/mL	50	s, p
Inhibin A	2.61 pg/mL	3 hrs+1 hr+15 min/RT; SH	0; 10-1000 pg/mL	50	s, p
Inhibin B Gen II	2.91 pg/mL	2+1 hr/RT; SH+30+10 min/RT; SH	0; 10-1000 pg/mL	50	s; p
Inhibin B Gen II Cal+CS	N.A.	N.A.	N.A.	N.A.	N.A.
Inhibin B Gen II RUO	2.91 pg/mL	2+1 hr/RT; SH+30+10 min/RT; SH	0; 10-1000 pg/mL	50	s; p
Inhibin B Gen II Cal+CS RUO	N.A.	N.A.	N.A.	N.A.	N.A.
Free Testosterone	<b>NEW</b> 0.097 pg/mL	2x1 hr/RT+10 min/RT	0; 0.1-75 pg/mL	50	s



## ALLERGY

Kit	Cat #	Wells	Enzyme	Wave-length
Histamine RUO	IM2015	12x8	AP	405
Histamine	IM2562	12x8	AP	405

Kit	Sensitivity	Incubation time	Calibrator range	Sample	
				Vol. [µL]	Form
Histamine RUO	0.52 nM	2 hrs/2-8°C+30 min/RT; SH*	0; 1-100 nM	100	p; u; b; c; l; ex:d
Histamine	0.52 nM	2 hrs/2-8°C+30 min/RT; SH*	0; 1-100 nM	100	p; u; b; c; l; ex:d

\*Non-shaking procedure available.



## DIABETES

Kit	Cat #	Wells	Enzyme	Wave-length
Anti-GAD	C15995	12x8	HRP	450+405
Anti-IA2	C15994	12x8	HRP	450+405

Kit	Sensitivity	Incubation time	Calibrator range	Sample	
				Vol. [µL]	Form
Anti-GAD	0.57 U/mL	2x1 hr/RT; SH+20 min/RT; SH+20 min/RT	0; 5-2000 U/mL	25	s
Anti-IA2	0.95 U/mL	(16-20 hrs+1)/2-8°C+20 min/RT; SH+20 min/RT	0; 7.5-4000 U/mL	50	s



## SPECIALTY

Kit	Cat #	Wells	Enzyme	Wave-length
25-OH Vitamin D total	B46328	12x8	HRP	450
1,25(OH) <sub>2</sub> Vitamin D	B68736	12x8	HRP	450
Histamine RUO	IM2369	12x8	AP	405
Serotonin RUO	IM1749	12x8	ACE	405
Chromogranin A	C24916	12x8	HRP	450
2-MET Plasma	C37563	2x(12x8)	HRP	450

Kit	Sensitivity (LoD*)	Incubation time	Calibrator range	Sample	
				Vol. [µL]	Form
25-OH Vitamin D total	2.81 ng/mL*	2 hrs+30+15 min/RT; SH	0; 10-100 ng/mL	50	s
1,25(OH) <sub>2</sub> Vitamin D	0.8 pg/mL	16-20 hrs/2-8°C+1 hr/2-8°C+15 min/RT	0; 3-180 pg/mL	50	s
Histamine RUO	0.06 µM	2x30 min/RT	0; 0.1-50 µM	50	ex:d (food)
Serotonin RUO	0.58 nM	30 min/RT+3 hrs+20 min/RT; SH	0; 1.6-200 nM	100	s; p; b; u; f; d
Chromogranin A	6.5 µg/L*	2x1 hr/RT; SH+30 min/RT; SH	0; 30-700 µg/L	50	s; p
2-MET Plasma	14.9/17.9 pg/mL*	15min/RT+15-20 hrs/2-8°C+30 min/RT+20-30 min/RT	0; 36-3600 pg/mL 0; 72-7200 pg/mL	200	ex: p

## AVAILABLE SEPARATELY

Component	Cat #	Intended for kits
AMH GEN II Assay Buffer (2x 26 mL)	B38561	A79765, A73818
TMB Chromogen Solution (11 mL)	DSL-10-9755-1	A79765, A73818, DSL-10-28100t-1, DSL-10-28100t-4, A81303, A81301
Wash Solution U (1 vial)	A54825	A79765, A73818, DSL-10-28100t-1, DSL-10-28100t-4, A81303, A81301
Inhibin B zero calibrator (5 mL)	C35899	A81303, A81301
Stop Solution (11 mL; H <sub>2</sub> SO <sub>4</sub> )	C24811	A79765, A73818, DSL-10-28100t-1, DSL-10-28100t-4, A81303, A81301
Stop Solution (4x 6mL; NaOH)	C03592	IM2369

ACE	Acetylcholinesterase	RUO	Research Use Only	SH	Shaking	ex:	Extraction of:
HRP	Horse Radish Peroxidase	N.A.	Not applicable	d	Solid biol. samples	s	Serum
AP	Alkaline Phosphatase	RT	Room temperature	f	Cerebrospinal fluid	p	Plasma
ELISA	Enzyme Linked Immunosorbent Assay	Cal	Calibrator	c	Cell culture supern.	b	Blood
LoD	Limit of Detection	CS	Control Sample	l	Other body liquids	u	Urine

Not all the products are available in all countries - please, ask your local distributor.

This product list is subject to change without notice.

All kits are IVD, CE marked, unless otherwise specified.



Labex Canada Inc.  
1091 Gorham St, Newmarket  
Ontario, L3Y8X7, Canada  
Phone: +1 (905) 392-7200  
Fax: +1 (905) 853-4548  
E-mail: [info@labexcanada.com](mailto:info@labexcanada.com)

Printed in the Czech Republic  
2020-10-26