

Primary Antibodies List - 2024

Product NO	Product name	Catalog number RTU	Catalog number Conc.	Antibody type	Clone	Cellular localization	Positive control	Antigen retrieving	Recommended dilutions	Certification(s)
0002	Actin, Muscle Specific (MSA)	LCR-0002-06	LCM-0002-1.0	Mouse monoclonal antibody	HHF35	Cytoplasm	Muscle/sarcoma	Citrate	50-100x	IVD CE
0004	Actin, Smooth Muscle (SMA)	LCR-0004-06	LCM-0004-1.0	Mouse monoclonal antibody	1A4	Cytoplasm	Smooth muscle/leiomyosarcoma	Citrate	100-200x	IVD CE
0007	ALK	LCR-0007-06	LCM-0007-1.0	Recombinant rabbit monoclonal antibody	rRM-A1	Cytoplasm and/or nucleus	Anaplastic lymphoma	Citrate	100-200x	IVD CE
0009	Alpha-1-Fetoprotein (AFP)	LCR-0009-06	LCM-0009-1.0	Rabbit polyclonal antibody	n/a	Cytoplasm	Fetal liver/hepatocellular carcinoma	EDTA	100-200x	IVD CE
0021	Androgen Receptor (AR)	LCR-0021-06	LCM-0021-1.0	Mouse monoclonal antibody	AR441	Nucleus	Prostate cancer	EDTA	100-200x	IVD CE
0024	Cyclin D1	LCR-0024-06	LCM-0024-1.0	Rabbit monoclonal antibody	SP4	Nucleus	Mantle cell lymphoma/breast cancer	Citrate	50-100x	IVD CE
0025	Bcl-2	LCR-0025-06	LCM-0025-1.0	Recombinant mouse monoclonal antibody	rMM-B1	Cellular membranes	Tonsillar/follicular lymphoma	EDTA	100-200x	Health Canada IVD CE
0060	Calcitonin (CT)	LCR-0060-06	LCM-0060-1.0	Rabbit monoclonal antibody	EP92	Cytoplasm	Thyroid carcinoma/medullary carcinoma	Citrate	50-100x	IVD CE
0061	Caldesmon	LCR-0061-06	LCM-0061-1.0	Recombinant rabbit monoclonal antibody	rRM-C1	Cytoplasm	Uterus	Citrate	100-200x	IVD CE
0062	Calponin	LCR-0062-06	LCM-0062-1.0	Recombinant rabbit monoclonal antibody	rRM-C2	Cytoplasm	Mammary ductal myoepithelial cells	EDTA	100-200x	IVD CE
0063	Calretinin (CR)	LCR-0062-06	LCM-0062-1.0	Rabbit monoclonal antibody	SP13	Cytoplasm	Mesothelioma	Citrate	100-200x	IVD CE
0069	CD1a	LCR-0069-06	LCM-0069-1.0	Recombinant mouse monoclonal antibody	rMM-C1	Cellular membranes	Thymus/tonsils	EDTA	50-100x	IVD CE
0075	CD3	LCR-0075-06	LCM-0075-1.0	Rabbit monoclonal antibody	SP41	Cellular membranes	Tonsils	EDTA	50-100x	Health Canada IVD CE
0080	CD5	LCR-0080-06	LCM-0080-1.0	Recombinant rabbit monoclonal antibody	rRM-C3	Cellular membranes	Tonsils	Citrate	50-100x	Health Canada IVD CE
0089	CD10	LCR-0089-06	LCM-0089-1.0	Recombinant mouse monoclonal antibody	rMM-C2	Cellular membranes	Tonsils/small intestine	EDTA	50-100x	Health Canada IVD CE

0099	CD19	LCR-0099-06	LCM-0099-1.0	Rabbit monoclonal antibody	SP110	Cellular membranes	Tonsils/lymph node	EDTA	100-200x	IVD CE
0100	CD20	LCR-0100-06	LCM-0100-1.0	Mouse monoclonal antibody	L26	Cellular membranes	Tonsils	Citrate	100-200x	Health Canada IVD CE
0101	CD21	LCR-0101-06	LCM-0101-1.0	Recombinant rabbit monoclonal antibody	rRM-C4	Cellular membranes	Tonsils	Citrate	50-100x	IVD CE
0104	CD23	LCR-0104-06	LCM-0104-1.0	Rabbit monoclonal antibody	EP75	Cellular membranes	Tonsils	Citrate	50-100x	IVD CE
0113	CD30/Ki-1	LCR-0113-06	LCM-0113-1.0	Recombinant rabbit monoclonal antibody	rRM-C5	Cellular membranes	Hodgkin's lymphoma	Citrate	50-100x	IVD CE
0114	CD31	LCR-0114-06	LCM-0114-1.0	Recombinant rabbit monoclonal antibody	rRM-C6	Cellular membranes	Tonsils/placenta	Citrate	50-100x	IVD CE
0117	CD34	LCR-0117-06	LCM-0117-1.0	Mouse monoclonal antibody	QBEnd/10	Cellular membranes	Tonsils	Citrate	50-100x	Health Canada IVD CE
0122	CD38	LCR-0122-06	LCM-0122-1.0	Rabbit monoclonal antibody	SN07-20	Cellular membranes	Tonsils	Citrate	100-200x	IVD CE
0126	CD43	LCR-0126-06	LCM-0126-1.0	Mouse monoclonal antibody	84-3C1	Cellular membranes	Tonsils/lymph node	EDTA	50-100x	IVD CE
0131	CD45/LCA	LCR-0131-06	LCM-0131-1.0	Mouse monoclonal antibody	Bra55	Cellular membranes	Tonsils	Citrate	100-200x	Health Canada IVD CE
0148	CD56	LCR-0148-06	LCM-0148-1.0	Mouse monoclonal antibody	56C04	Cellular membranes	Neuroblastoma/ retinoblastoma	Citrate	100-200x	IVD CE
0160	CD68	LCR-0160-06	LCM-0160-1.0	Recombinant mouse monoclonal antibody	rMM-C3	Cytoplasm	Tonsils/lymph node/spleen	Citrate	100-200x	IVD CE
0165	CD79a	LCR-0165-06	LCM-0165-1.0	Mouse monoclonal antibody	HM47/ArMM-C3	Cellular membranes	Tonsils/lymph node	Citrate	100-200x	Health Canada IVD CE
0175	CD99	LCR-0175-06	LCM-0175-1.0	Mouse monoclonal antibody	HO36-1.1	Cellular membranes	Ewing's sarcoma/pancreas	EDTA	25-50x	IVD CE
0184	CD117/c-kit	LCR-0184-06	LCM-0184-1.0	Rabbit monoclonal antibody	EP10	Cellular membranes/cytoplasm	Gastrointestinal stromal tumor/tonsils/skin	Citrate	50-100x	Health Canada IVD CE
0194	Carcinoembryonic Antigen (CEA)	LCR-0194-06	LCM-0194-1.0	Mouse monoclonal antibody	COL-1	Cytoplasm	Colon cancer	Citrate	50-100x	IVD CE
0196	c-erbB-2/HER2	LCR-0196-06	LCM-0196-1.0	Rabbit monoclonal antibody	EP3	Cellular membranes	Breast cancer	Citrate	50-100x	Health Canada IVD CE

0201	HCG-beta	LCR-0201-06	LCM-0201-1.0	Mouse monoclonal antibody	HCGb/54+ HCGb/459	Cytoplasm	Placental/choriocarcinoma	Citrate	50-100x	IVD CE
0202	Chromogranin A (CgA)	LCR-0202-06	LCM-0202-1.0	Recombinant mouse monoclonal antibody	rMM-C4	Cytoplasm	Pancreas	Citrate	100-200x	Health Canada IVD CE
0231	Desmin	LCR-0231-06	LCM-0231-1.0	Recombinant rabbit monoclonal antibody	rRM-D1	Cytoplasm	Muscle/sarcoma	Citrate	50-100x	IVD CE
0236	Epithelial Membrane Antigen (EMA)	LCR-0236-06	LCM-0236-1.0	Mouse monoclonal antibody	GP1.4	Cytoplasm and cellular membranes	Breast/colon cancer	Citrate	100-200x	IVD CE
0238	EGFR	LCR-0238-06	LCM-0238-1.0	Recombinant rabbit monoclonal antibody	rRM-E1	Cellular membranes	Squamous cell carcinoma/skin/placenta	EDTA	25-50x	Health Canada IVD CE
0241	Estrogen Receptor (ER)	LCR-0241-06	LCM-0241-1.0	Rabbit monoclonal antibody	SP1	Nucleus	Breast cancer	EDTA	50-100x	Health Canada IVD CE
0255	Follicle Stimulating Hormone (FSH)	LCR-0255-06	LCM-0255-1.0	Mouse monoclonal antibody	FSH03	Cytoplasm	Anterior pituitary	Citrate	50-100x	IVD CE
0257	Galectin-3	LCR-0257-06	LCM-0257-1.0	Mouse monoclonal antibody	B2C10	Cytoplasm	Anaplastic large cell lymphoma/papillary thyroid carcinoma	Citrate	50-100x	IVD CE
0262	GFAP	LCR-0262-06	LCM-0262-1.0	Recombinant mouse monoclonal antibody	rMM-G1	Cytoplasm	Brain/astrocytoma	Citrate	50-100x	IVD CE
0279	Human Growth Hormone (hGH)	LCR-0279-06	LCM-0279-1.0	Mouse monoclonal antibody	GH01 #2 +GH02 #7	Cytoplasm	Recombinant human growth hormone /human anterior pituitary	Citrate	50-100x	IVD CE
0282	HBcAg	LCR-0282-06	LCM-0282-1.0	Rabbit polyclonal antibody	n/a	Nucleus/perinuclear	Hepatitis B virus infected liver cells	Citrate	100-200x	IVD CE
0283	HBsAg	LCR-0283-06	LCM-0283-1.0	Mouse monoclonal antibody	A10F1	Cytoplasm	Hepatitis B virus infected liver cells	EDTA	100-200x	IVD CE
0291	Helicobacter phlori (HP)	LCR-0291-06	LCM-0291-1.0	Rabbit polyclonal antibody	n/a	Bacteria on the surface of the gastric mucosa	Gastric mucosa infected with Helicobacter pylori	EDTA	100-200x	Health Canada IVD CE
0332	Cytokeratin7 (CK7)	LCR-0332-06	LCM-0332-1.0	Recombinant mouse monoclonal antibody	rMM-C5	Cytoplasm	Ovarian cancer	Citrate	100-200x	Health Canada IVD CE
0333	Cytokeratin8 (CK8)	LCR-0333-06	LCM-0333-1.0	Mouse monoclonal antibody	TS1	Cytoplasm	Skin/colon/ lung/breast cancer	Citrate	100-200x	IVD CE
0334	Cytokeratin8/18 (CK8/18)	LCR-0334-06	LCM-0334-1.0	Mouse monoclonal antibody	K8.8+ DC10	Cytoplasm	Skin/Breast cancer	EDTA	50-100x	IVD CE
0341	Cytokeratin18 (CK18)	LCR-0341-06	LCM-0341-1.0	Mouse monoclonal antibody	DC10	Cytoplasm	Skin/breast cancer/lung cancer	Citrate	100-200x	IVD CE

0342	Cytokeratin19 (CK19)	LCR-0342-06	LCM-0342-1.0	Recombinant rabbit monoclonal antibody	rRM-C9	Ccytoplasm	Skin/colon cancer	EDTA	100-200x	IVD CE
0343	Cytokeratin20 (CK20)	LCR-0343-06	LCM-0343-1.0	Recombinant rabbit monoclonal antibody	rRM-C10	Cytoplasm	Colon cancer	Citrate	100-200x	Health Canada IVD CE
0346	Cytokeratin, HMW (HCK)	LCR-0346-06	LCM-0346-1.0	Mouse monoclonal antibody	34 β E12	Cytoplasm	Squamous cell carcinoma	Citrate	50-100x	IVD CE
0347	Cytokeratin, LMW (LCK)	LCR-0346-06	LCM-0346-1.0	Mouse monoclonal antibody	AE1	Cytoplasm	Skin/squamous cell carcinoma	EDTA	50-100x	IVD CE
0349	Cytokeratin, Pan (CK Pan)	LCR-0347-06	LCM-0347-1.0	Recombinant mouse monoclonal antibody	rMM-C7	Cytoplasm	Skin/Lung cancer	EDTA	50-100x	Health Canada IVD CE
0350	Ki-67	LCR-0350-06	LCM-0350-1.0	Recombinant mouse monoclonal antibody	rMM-K1	Nucleus	Tonsils/breast cancer/spleen	EDTA	50-100x	Health Canada IVD CE
0373	MART-1/Melan-A	LCR-0373-06	LCM-0373-1.0	Recombinant mouse monoclonal antibody	rMM-M1	Cytoplasm	Melanoma	Citrate	50-100x	Health Canada IVD CE
0379	Melanoma (HMB45)	LCR-0379-06	LCM-0379-1.0	Mouse monoclonal antibody	HMB45	Cytoplasm	Melanoma	Citrate	50-100x	Health Canada IVD CE
0380	Mesothelial Cell (MC)	LCR-0379-06	LCM-0379-1.0	Recombinant mouse monoclonal antibody	rMM-M2	Cytoplasm	Mesothelioma	Citrate	50-100x	IVD CE
0396	MUC-1	LCR-0396-06	LCM-0396-1.0	Mouse monoclonal antibody	MUC1/845	Cytoplasm and cellular membranes	Breast cancer	Citrate	100-200x	Health Canada IVD CE
0397	MUC-2	LCR-0397-06	LCM-0397-1.0	Recombinant rabbit monoclonal antibody	rRM-M1	Cytoplasm and cellular membranes	Small intestine	Citrate	100-200x	IVD CE
0399	Mucin 5AC (MUC5AC)	LCR-0399-06	LCM-0399-1.0	Mouse monoclonal antibody	ZM148	Cytoplasm	Stomach	Citrate	50-100x	IVD CE
0405	Myeloperoxidase (MPO)	LCR-0405-06	LCM-0405-1.0	Rabbit monoclonal antibody	EP151	Cytoplasm	Tonsils/spleen	Citrate	100-200x	IVD CE
0409	MyoD1	LCR-0409-06	LCM-0409-1.0	Rabbit monoclonal antibody	EP212	Nucleus	Rhabdomyosarcoma	EDTA	50-100x	IVD CE
0410	Myogenin	LCR-0410-06	LCM-0410-1.0	Rabbit monoclonal antibody	F5D	Nucleus	Rhabdomyosarcoma/ Wilms sarcoma	Citrate	50-100x	IVD CE
0415	Neurofilament (NF)	LCR-0415-06	LCM-0415-1.0	Recombinant rabbit monoclonal antibody	rRM-N1	Cytoplasm	Brain	Citrate	50-100x	IVD CE
0416	Neuron Specific Enolase (NSE)	LCR-0416-06	LCM-0416-1.0	Mouse monoclonal antibody	ENO2/1375	Cytoplasm	Intestine/Pancreas	Citrate	100-200x	IVD CE

0421	nm23	LCR-0421-06	LCM-0421-1.0	Mouse monoclonal antibody	12A1	Cytoplasm and nucleus	Breast cancer	EDTA	100-200x	IVD CE
0425	p16	LCR-0425-06	LCM-0425-1.0	Recombinant mouse monoclonal antibody	rRM-P1	Cytoplasm/nucleus	Cervical squamous cell carcinoma	Citrate	50-100x	Health Canada IVD CE
0430	p53	LCR-0430-06	LCM-0430-1.0	Recombinant mouse monoclonal antibody	rRM-P2	Nucleus	Colon/Breast cancer	Citrate	50-100x	IVD CE
0442	Placental Alk Phosphatase (PLAP)	LCR-0442-06	LCM-0442-1.0	Recombinant rabbit monoclonal antibody	rRM-P1	Cellular membranes	Placenta/seminoma	Citrate	100-200x	Health Canada IVD CE
0443	Human Placental Lactogen(hPL)	LCR-0443-06	LCM-0443-1.0	Recombinant rabbit monoclonal antibody	rRM-H1	Cytoplasm	Placenta	Citrate	50-100x	IVD CE
0448	Progesterone Receptor (PR)	LCR-0448-06	LCM-0448-1.0	Rabbit monoclonal antibody	SP2	Nucleus	Breast cancer	EDTA	100-200x	Health Canada IVD CE
0453	Prostate Specific Antigen (PSA)	LCR-0453-06	LCM-0453-1.0	Rabbit monoclonal antibody	EP109	Cytoplasm	Normal prostate/prostate cancer	Citrate	100-200x	Health Canada IVD CE
0455	Prostate Specific Acid Phosphatase (PSAP)	LCR-0455-06	LCM-0455-1.0	Mouse monoclonal antibody	PASE/4LJ	Cytoplasm	Normal prostate/prostate cancer	Citrate	100-200x	IVD CE
0466	S-100 Protein	LCR-0466-06	LCM-0466-1.0	Recombinant mouse monoclonal antibody	rMM-S1	Cytoplasm	Melanoma/neurilemmoma (Schwannoma)	Citrate	100-200x	Health Canada IVD CE
0482	Synaptophysin (Syn)	LCR-0482-06	LCM-0482-1.0	Rabbit monoclonal antibody	EP158	Cytoplasm	Pancreas/pheochromocytoma	Citrate	100-200x	IVD CE
0489	Terminal Deoxyribonucleotidyl Transferase (TdT)	LCR-0489-06	LCM-0489-1.0	Rabbit monoclonal antibody	EP266	Nucleus	Thymus	Citrate	25-50x	IVD CE
0495	Thyroglobulin (TG)	LCR-0495-06	LCM-0495-1.0	Mouse monoclonal antibody	TGB04+ TGB05	Cytoplasm	Thyroid gland	Citrate	100-200x	IVD CE
0497	Thyroid Stimulating Hormone (TSH)	LCR-0497-06	LCM-0497-1.0	Mouse monoclonal antibody	TSH01+ TSH02	Cytoplasm	Anterior pituitary	Citrate	50-100x	IVD CE
0503	Topoisomerase II alpha (TOP2A)	LCR-0503-06	LCM-0503-1.0	Recombinant mouse monoclonal antibody	rRM-T1	Nucleus	Tonsils	Citrate	100-200x	IVD CE
0504	p504s/AMACR	LCR-0504-06	LCM-0504-1.0	Mouse monoclonal antibody	13H4	Cytoplasm	Prostate cancer	EDTA	100-200x	IVD CE
0514	Thyroid TranscriptionFactor-1 (TTF-1)	LCR-0514-06	LCM-0514-1.0	Mouse monoclonal antibody	8G7G3/1	Nucleus	Thyroid/Lung	EDTA	50-100x	Health Canada IVD CE
0523	Vascular Endothelial Growth Factor (VEGF)	LCR-0523-06	LCM-0523-1.0	Mouse monoclonal antibody	VG1	Cytoplasm, cell surface, extracellular matrix	Angiosarcoma	EDTA	25-50x	IVD CE

0524	Vimentin	LCR-0524-06	LCM-0524-1.0	Rabbit monoclonal antibody	SP20	Cytoplasm	Sarcoma/melanoma	Citrate	100-200x	IVD CE
0526	Wilm's Tumor (WT1)	LCR-0526-06	LCM-0526-1.0	Mouse monoclonal antibody	WT1/857	Nucleus	Wilms' tumor	EDTA	25-50x	IVD CE
0536	E-Cadherin	LCR-0536-06	LCM-0536-1.0	Recombinant mouse monoclonal antibody	rMM-E1	Cellular membranes	Breast/colon cancer	EDTA	50-100x	Health Canada IVD CE
0543	CD138	LCR-0543-06	LCM-0543-1.0	Mouse monoclonal antibody	MCA2459	Cytoplasm & cellular membraness	Tonsils	EDTA	100-200x	IVD CE
0545	Catenin (beta)	LCR-0545-06	LCM-0545-1.0	Mouse monoclonal antibody	CTNNB1/1508	Cellular membranes	Breast cancer	Citrate	50-100x	Health Canada IVD CE
0568	Bcl-6	LCR-0568-06	LCM-0568-1.0	Recombinant mouse monoclonal antibody	rRB-B2	Nucleus	Tonsils	EDTA	25-50x	IVD CE
0589	Villin	LCR-0589-06	LCM-0589-1.0	Mouse monoclonal antibody	CWWB1	Cytoplasm cellular membraness	Colon/rectum	Citrate	50-100x	IVD CE
0640	MSH2	LCR-0640-06	LCM-0640-1.0	Rabbit monoclonal antibody	RED2	Nucleus	Tonsils	Citrate	100-200x	Health Canada IVD CE
0654	p63	LCR-0654-06	LCM-0654-1.0	Recombinant rabbit monoclonal antibody	rRM-P2	Nucleus	Tonsils	Citrate	50-100x	Health Canada IVD CE
0658	Perforin	LCR-0658-06	LCM-0658-1.0	Mouse monoclonal antibody	E5	Cytoplasm	Spleen	Citrate	100-200x	IVD CE
0665	PTEN	LCR-0665-06	LCM-0665-1.0	Rabbit polyclonal antibody	n/a	Cytoplasm	Tonsils	Citrate	50-100x	Health Canada IVD CE
0669	Surfactant Protein B (SP-B)	LCR-0669-06	LCM-0669-1.0	Mouse monoclonal antibody	F-2	Cytoplasm	Pulmonary adenocarcinoma	Citrate	50-100x	IVD CE
0688	Survivin	LCR-0688-06	LCM-0688-1.0	Rabbit monoclonal antibody: EP119	EP119	Nucleus/cytoplasm	Stomach/normal pancreas/prostate cancer	Citrate	50-100x	IVD CE
0690	MLH1	LCR-0690-06	LCM-0690-1.0	Mouse mab	4C9C7	Nucleus	Colon cancer	EDTA	100-200x	Health Canada IVD CE
0713	CDX-2	LCR-0713-06	LCM-0713-1.0	Rabbit monoclonal antibody	EPR2764Y	Nucleus	Colon cancer	Citrate	50-100x	Health Canada IVD CE
0715	COX-2	LCR-0715-06	LCM-0715-1.0	Rabbit monoclonal antibody	SP21	Cytoplasm and cellular membranes	Colorectal adenocarcinoma	Citrate	25-50x	IVD CE
0740	MUM1	LCR-0740-06	LCM-0740-1.0	Rabbit monoclonal antibody	EP190	Nucleus	Tonsils/spleen	EDTA	50-100x	IVD CE

0741	Thyroid Peroxidase (TPO)	LCR-0741-06	LCM-0741-1.0	Rabbit monoclonal antibody	EP159	Cytoplasm	Thyroid gland	EDTA	100-200x	IVD CE
0744	Surfactant Protein-A (SP-A)	LCR-0744-06	LCM-0744-1.0	Mouse monoclonal antibody	6F10	Cytoplasm	Lung cancer (non-mucinous bronchioloalveolar carcinoma)	Citrate	25-50x	IVD CE
0745	Cytokeratin5/6 (CK5/6)	LCR-0745-06	LCM-0745-1.0	Recombinant mouse monoclonal antibody	rMM-C8	Cytoplasm	Prostate, tonsils, skin tissue	Citrate	50-100x	Health Canada IVD CE
0749	Inhibin, Alpha	LCR-0749-06	LCM-0749-1.0	Rabbit monoclonal antibody	EP378	Cytoplasm	Adrenal cortex, placenta, testis	Citrate	25-50x	IVD CE
0754	D2-40	LCR-0754-06	LCM-0754-1.0	Mouse monoclonal antibody	D2-40	Cytoplasm	Breast cancer, mesothelioma and colon cancer	EDTA	100-200x	IVD CE
0758	p120	LCR-0758-06	LCM-0758-1.0	Mouse monoclonal antibody	CTTNB1/1508	Cellular membranes and cytoplasm	Ductal invasive carcinoma of the breast	Citrate	50-100x	IVD CE
0760	DOG-1	LCR-0760-06	LCM-0760-1.0	Mouse monoclonal antibody	DOG1.1	Cytoplasm & cellular membraness	Gastrointestinal stromal tumors	Citrate	100-200x	IVD CE
0777	Napsin A	LCR-0777-06	LCM-0777-1.0	Mouse monoclonal antibody	ZM11	Cytoplasm	Pulmonary adenocarcinoma	Citrate	100-200x	Health Canada IVD CE
0778	p40	LCR-0778-06	LCM-0778-1.0	Recombinant rabbit monoclonal antibody	rRM-P3	Nucleus	Pulmonary squamous cell carcinoma	EDTA	100-200x	Health Canada IVD CE
0779	MSH6	LCR-0779-06	LCM-0779-1.0	Rabbit monoclonal antibody	EP49	Nucleus	Colon cancer	Citrate	50-100x	Health Canada IVD CE
0780	PMS2	LCR-0780-06	LCM-0780-1.0	Rabbit monoclonal antibody: EP51	EP51	Nucleus	Colon	EDTA	50-100x	Health Canada IVD CE
0781	GATA3	LCR-0781-06	LCM-0781-1.0	Recombinant mouse monoclonal antibody	rMM-G1	Nucleus	Breast, bladder cancer tissue	Citrate	50-100x	IVD CE
0782	Mammaglobin	LCR-0782-06	LCM-0782-1.0	Mouse monoclonal antibody	304-1A5	Cytoplasm	Breast cancer tissue	Citrate	50-100x	IVD CE
0783	Pax-8	LCR-0783-06	LCM-0783-1.0	Mouse monoclonal antibody: 6G5	6G5	Nucleus	Renal cell carcinoma	EDTA	100-200x	Health Canada IVD CE
0785	PD-1	LCR-0785-06	LCM-0785-1.0	Mouse monoclonal antibody	2E5	Cellular membranes	Angioimmunoblastic T-cell lymphoma/tonsils	EDTA	50-100x	IVD CE
0786	PD-L1	LCR-0786-06	LCM-0786-1.0	Rabbit monoclonal antibody: 1H3	1H3	Cellular membranes	Placenta or lung's adenocarcinoma	EDTA	100-200x	IVD CE

- NOTES:**
- **Only products with MDL licences, registered by Health Canada (marked "Health Canada" in the "Certifications" column) can be sold in Canada.**
 - **Not all the products are available in all countries, please ask your local distributor.**